

Double action IMT/PEP DEVICE

Ordering code: BHC-IMTPEP-1

Ask us
about
supply!



TYPE 1 IMT MODE

Diaphragm expansion

Strengthens
the diaphragm



TYPE 2 PEP MODE

Diaphragm contraction

Strengthens
exhalation



TYPE 3 DOUBLE ACTION

IMT + PEP

Simultaneous contraction
and expansion

Use the IMT/PEP Device to:

- Clear the throat
- Draw out mucus
- Develop stronger breathing

IMT pressure range - 10-40 cmH2O
PEP pressure range- 5-20 cmH2O



Increase Lung Capacity



Strengthen Respiratory Muscles



Increase Cardiac Output

Trustworthy

We currently supply our products to local hospitals and health centers for breathing rehabilitation.

Easy to Clean

To clean the device, remove the mouthpiece and the PEP body and wash everything using warm running water. Dry everything well so the spring inside the body will not rust.

Double Action

The IMT/PEP device is a three-mode breathing exerciser effective in breathing exercises and rehabilitation. It enables easier mucus removal through exhalation and expands lung capacity through inhalation.

Performs each exercise separately or simultaneously.